

Collection Device



This list may not contain all physician requests and acceptable substitute collection supplies may be provided when devices shown are not available. Please refer to the online Aspirus Reference Laboratory Test Directory (testmenu.com/Aspirus) for the most comprehensive and up to date specimen collection information.

Test

- Culture, Aerobic
- Culture, Anaerobic with Aerobe
- Culture, Yeast
- Culture, Eye
- Culture, Gonococcus
- Culture, Ear
- Staph aureus Culture
- Culture, Vancomycin-Resistant Enterococcus
- Culture, Methicillin Resistant Staphylococcus aureus
- Group B Strep NAT, Amplified
- Culture, Fungus
- Fungal Smear
- Staph aureus & MRSA Detection by NAT
- Wet Prep, Genital

Collection Device - ESwab Regular



Test

- Bordetella pertussis/parapertussis, NAT Amplified

Collection Device - SWAB-PAK



Test

- Culture, Nasopharyngeal
- Pediatric patient cultures

Collection Device - ESwab Flex Mini



Test

- Vaginal Panel NAT, Amplified

Collection Device - BD Molecular Swab Kit



Test

- Rapid Strep Group A with Reflex to Culture


Collection Device - Dry Swab - for Rapid Strep ESwab Regular - for Culture



Test

- Viral Culture
- Culture, Herpes Simplex Virus
- Culture, Viral, Non-Respiratory
- Culture, Chlamydia trachomatis
- Culture, Mycoplasma & Ureaplasma

Collection Device - ARUP Transport Media & Swab Kit



Test

- Chlamydia trachomatis NAT, Amplified
- Neisseria gonorrhoeae NAT, Amplified
- Trichomonas vaginalis NAT, Amplified
- Eye/Oral/Rectal: Miscellaneous Lab Test, Modifier 0060241

Collection Device - Aptima Collection Kits

Source	Collection Kit	
<ul style="list-style-type: none">• Provider or Patient collected vaginal swab• Provider collected rectal swab (not for T. vag NAT)• Provider collected pharyngeal swab (not for T. vag NAT)	Aptima Multitest Swab Collection Kit	
	Aptima Unisex Collection Kit (Use the blue shaft swab for sample collection.)	
	Aptima Urine Collection Kit <ul style="list-style-type: none">• Transfer approximately 2ml of first-catch urine into Aptima urine collection tube.• Urine volume must be between the black indicators on the tube.• Minimum volume = 2 ml	
Liquid Cytology	Aptima Specimen Transfer Kit <ul style="list-style-type: none">• 1 ml of a well-mixed ThinPrep or SurePath liquid cytology specimen transferred to an Aptima Specimen Transfer Tube.	

Test

- Influenza A/B, NAT, Amplified
- Influenza A & B Rapid Antigen
- Influenza A/B, RSV, NAT, Amplified
- SARS-CoV2, NAT, Amplified
- SARS-CoV-2 & Influenza A/B NAT, Amplified
- SARS-CoV-2, Influenza A/B, & RSV NAT, Amplified
- RSVE
- RSV, NAT, Amplified

Collection Device - Nasal or NP Swab UTM (Universal Transport Media)

